


CLONOSEQ MRD TESTING IN ADULT ACUTE LYMPHOBLASTIC LEUKEMIA (ALL): A POTENTIAL CLINICAL PATHWAY

Based on available evidence, the following clinical strategy reflects a potential approach for integrating clonoSEQ® MRD testing into the management of patients with adult ALL:

Key Takeaways

- Clinical practice guidelines suggest performing a Clonality (ID) assessment at time of diagnosis to establish tumor-specific DNA sequences to track; use a fresh or archived bone marrow or blood sample collected prior to initiation of treatment.¹
- Evaluation of MRD status should be considered after induction, after consolidation, and at further time points during maintenance and follow-up.^{2,3,4,5,6}
- In general, MRD positivity at the end of induction/consolidation predicts high relapse rates and should prompt an evaluation for allogeneic hematopoietic stem cell transplantation (HSCT).^{1,4,7}
- Prospective MRD monitoring is recommended for MRD-negative patients in 3 to 6-month intervals during maintenance therapy and follow-up.^{1,3,6,8}



Supporting Data, Guidelines and References

The study design above reflects guideline-recommended MRD timepoints and incorporates assessment timelines from the 06/99 and 07/03 trials of the German Multicenter Study Group for Adult Acute Lymphoblastic Leukemia (GMALL).³ In addition to guidance for testing during treatment, data generated from these trials and related studies support the use of serial MRD monitoring in adult ALL post-transplant and during surveillance.

MRD Assessment Timepoints

- Per NCCN guidelines, it is recommended that the initial MRD measurement be performed on completion of induction therapy; additional time points for MRD evaluation should be guided by the treatment protocol or regimen used.¹

- In line with GMALL 06/99 and 07/03, consider MRD evaluation after induction (day 71), after consolidation (week 16), and in 3-to-6-month intervals during maintenance therapy and follow-up.^{2,3,4}
 - *Note: these testing recommendations are presented within the context of the specific 6-drug regimen utilized in GMALL 06/99 and 07/03. Precise timepoints for MRD testing when using other regimens will require modification accordingly.*
- Consider HSCT for those patients deemed MRD-positive after consolidation.^{4,7}

Summary of supporting data:

- Studies in adults with ALL have shown a strong correlation between MRD positivity and risk for relapse.^{1,2,5,6,9,10,11}
- Existing evidence suggests high value of MRD monitoring for identifying subgroups of patients who may benefit from further intensified therapies or alternative treatments.
 - Data supports unique risk stratification via MRD assessment after induction, then again after consolidation.^{2,3,4,5}
 - Post-consolidation MRD status is a risk factor for poorer outcomes in adults with ALL and may stratify high-risk patients who might benefit from HSCT.^{4,7}
- Serial monitoring of MRD positivity may be useful in patients with molecular relapse and low-level disease after induction and consolidation, and patients with persistent MRD negativity have been shown to have a favorable OS rate after HSCT.^{1,6,8,12}
- MRD negativity via NGS (10^{-6}) has been shown to predict relapse in adult patients receiving allo-HSCT, and may confer additional prognostic information compared to MRD status by less sensitive measures (e.g., MFC) at identical time points.^{12,13,14}
- Studies have shown a high correlation between MRD results assessed in the bone marrow and peripheral blood in the context of adult ALL; thus, blood* may be a potential sample source for MRD assessment under certain circumstances.^{12,15,16}

Summary of guidelines:

- Guidelines recommend MRD following induction therapy, with the optimal sample for MRD assessment being the first pull or early pull of the bone marrow aspirate.¹
- For patients with unavailable MRD status at start of consolidation, consider retesting for MRD at first available opportunity.¹
- When possible, therapy aimed at eliminating MRD prior to allogeneic HSCT is preferred.¹

1. NCCN Clinical Practice Guidelines in Oncology: Acute Lymphoblastic Leukemia. Version 2.2021.
2. Brüggemann M, et al. *Blood*. 2006;107(3):1116-1123.
3. Raff T, et al. *Blood*. 2007;109(3):910-915.
4. Gökbuget N, et al. *Blood*. 2012;120(9):1868-1876.
5. Vidrales MB, et al. *Blood*. 2003;101(12):4695-4700.
6. Patel B, et al. *Br J Haematol*. 2010;148(1):80-89.
7. Bassan R, et al. *Blood*. 2009;113(18):4153-4162.
8. Pemmaraju N, et al. *Am J Hematol*. 2017;92(3):279-285.
9. Mortuza FY, et al. *J Clin Oncol*. 2002;20(4):1094-1104.
10. Holowiecki J, et al. *Br J Haematol*. 2008;142(2):227-237.
11. Ravandi F, et al. *Br J Haematol*. 2016;172(3):392-400.
12. Logan AC, et al. *Biol Blood Marrow Transplant*. 2014;20(9):1307-1313.
13. Mannis GN, et al. *Biol Blood Marrow Transplant*. 2016;22(6):1030-1036.
14. Short NJ et al. ASH 2020; Abstract 583.
15. Sala Torra O, et al. *Biol Blood Marrow Transplant*. 2017;23(4):691-696.
16. Muffly LS, et al. *Blood*. 2020;136(Suppl 1):22-23.

*Blood-based clonoSEQ testing in ALL is available as a CLIA-validated LDT and has not been cleared or approved by the FDA.

clonoSEQ® is available as an FDA-cleared *in vitro* diagnostic (IVD) test service provided by Adaptive Biotechnologies to detect minimal residual disease (MRD) in bone marrow from patients with multiple myeloma or B-cell acute lymphoblastic leukemia (B-ALL) and blood or bone marrow from patients with chronic lymphocytic leukemia (CLL). Additionally, clonoSEQ is available for use in other lymphoid cancers and specimen types as a CLIA validated laboratory developed test (LDT). For important information about the FDA-cleared uses of clonoSEQ including test limitations, please visit clonoSEQ.com/technical-summary.